

InBody



InBody270s

Body Composition Analyzer for All

Body composition analyzer designed for everyday use by everyone

New Standard of BIA

Redefine the depth of health assessment with the Phase Angle

Streamlined Portability


Compact and foldable design to ensure effortless mobility

InBody Technology

InBody uses Bioelectrical Impedance Analysis (BIA) technology to measure human body composition. Impedance is the resistance of the human body generated when a micro alternating current flows through the human body. The human body is made of water that conducts electricity well, and the resistance varies depending on the amount of water. BIA is a technology that quantitatively measures body water through impedance that occurs when an electric current flows through the human body. InBody provides diverse information on body composition based on the measured body water.

Direct Segmental Measurement-BIA

The human body exhibits varying lengths and cross-sectional areas for each body segments. Arms and legs, characterized by narrow cross-sectional areas and length, exhibit higher impedance values and lower muscle mass. Conversely, the trunk, with its broader cross-sectional area, yields lower impedance values and higher muscle mass. Even the slightest change in trunk impedance can significantly influence the total muscle mass. Therefore, it is essential to separately measure trunk impedance for precise total muscle mass assessment. InBody conducts separate measurements for arms, legs, and the trunk, ensuring the utmost accuracy in the analysis.




8-Point Tactile Electrodes utilizing Thumb Electrodes

Using the structural features of the human body, InBody pioneered '8-Point Tactile electrode with Thumb Electrodes'. This ensures InBody measurements start at the same location on the wrists and ankles, guaranteeing reliable and reproducible results.




Simultaneous Multi-Frequency Impedance Measurement

InBody introduced a technology in body composition analyzers to transmit multiple frequencies at once, obtaining specific impedance data for each for the first time. This reduces measurement time and error, leading to more accurate body water and fluid balance measurements.




No Estimations or Empirical Equations on Measured Values

InBody does not rely on empirical estimations based on age, gender, and more to ensure the accuracy of the measured data. In the past, empirical estimations were applied to the equations to ensure accuracy due to technological limitations. However, this resulted in lower accuracy when the measured population group changes. InBody overcame these limitations with technological developments such as direct segmental measurement-BIA to measure and analyze accurate body composition without applying empirical estimation. Therefore, InBody devices can provide data regardless of population and can reflect changes in the body with higher sensitivity.



Over 98% Correlation to DEXA on Accuracy

InBody precisely detects changes in body composition using impedance alone, showing a correlation over 0.98 with the gold-standard DEXA device.



Ryan T Hurt et al., The Comparison of SMF-BIA and DEXA for Estimating Fat Free Mass and Percentage Body Fat in an Ambulatory Population, *J Parenter Enteral Nutr.* 2021 Aug;45(6):1231-1238

InBody270S Highlights

Quick Measurement

Experience InBody's quick and precise measurement in only 30 seconds, with ensured reliability. Users can access their health data instantly, facilitating prompt consultations.

Compact and Portable Design

The foldable structure and compact size ensure easy transportation and spatial efficiency. The optional InBody270 carrying bag enhances portability. Its compact size maximizes space utilization without compromising on functionality.

Results in Your Hands

Transfer your body composition data directly to your mobile phone using our QR code feature. Stay connected and track your progress anytime, anywhere with instant access to your detailed results.



Comprehensive Parameters for All

Nutrition Assessment

In-depth Body Composition Analysis

Nutrition Analysis offers comprehensive evaluations of protein, mineral levels, and body fat to ensure the body's nutritional needs are met for optimal health. By providing insights into protein consumption, mineral status, and percent body fat, it empowers users to make informed dietary choices. With this in-depth analysis, achieving a balanced diet and monitoring progress toward health goals becomes straightforward and manageable.

Cellular Integrity Check

Phase Angle

The human body comprises 36 trillion cells, and understanding cell health is crucial for overall well-being. The Phase Angle is a key parameter in assessing cell health and overall physiological status. It reflects the relationship between resistance in total body water and reactance in cell membrane. A higher Phase Angle indicates better cell membrane integrity, and well-balanced fluid, suggesting healthier cells. Last but not least, with the addition of Whole Body Phase Angle History, users can intuitively track and monitor their health trends over time.

Sarcopenia Assessment

SMI(Skeletal Muscle Mass Index)

Sarcopenia, assigned the diagnosis code M62.84 by WHO, is acknowledged as a disease rather than just a natural phenomenon. It can be easily assessed and evaluated using the Skeletal Muscle Mass Index (SMI)* and Hand Grip Strength**, allowing for comprehensive evaluation and personalized consultations.

*Skeletal Muscle Mass Index (SMI) calculated by taking the sum of the appendicular muscle mass (in kilograms) and dividing it by the square of the person's height (in meters).

**Hand Grip Strength is available with connections to the InBody Handgrip Dynamometer (IB-HGS, optional).



InBody Result Sheet

Provides reference parameters to thoroughly evaluate patients' conditions across various medical practices.

InBody

[InBody270S]

7

Customized Logo

www.customized.com

ID	Height	Age	Gender	Test Date & Time
Jane Doe	156.9cm	51	Female	05.30.2024 11:16

1 Body Composition Analysis

Total amount of water in body	Total Body Water	(L)	28.1 (27.0 ~ 33.0)
For building muscles	Protein	(kg)	7.4 (7.2 ~ 8.8)
For strengthening bones	Minerals	(kg)	2.76 (2.49 ~ 3.05)
For storing excess energy	Body Fat Mass	(kg)	20.8 (10.6 ~ 16.9)
Sum of the above	Weight	(kg)	59.1 (45.0 ~ 60.8)


2 Muscle-Fat Analysis

	Under	Normal	Over
Weight (kg)	55 70 85 100 115 130 145 160 175 190 205 %	59.1	
SMM (kg) <small>Skeletal Muscle Mass</small>	70 80 90 100 110 120 130 140 150 160 170 %	20.4	
Body Fat Mass (kg)	40 60 80 100 160 220 280 340 400 460 520 %	20.8	


3 Obesity Analysis

	Under	Normal	Over
BMI (kg/m ²) <small>Body Mass Index</small>	10.0 15.0 18.5 21.5 25.0 30.0 35.0 40.0 45.0 50.0 55.0	24.0	
PBF (%) <small>Percent Body Fat</small>	8.0 13.0 18.0 23.0 28.0 33.0 38.0 43.0 48.0 53.0 58.0	35.2	

4 Segmental Lean Analysis



5 Segmental Fat Analysis



* Segmental fat is estimated.

6 Body Composition History

Weight	65.3	63.9	62.4	61.8	62.3	60.9	60.5	59.1
SMM (kg) <small>Skeletal Muscle Mass</small>	20.1	20.0	19.7	19.7	19.8	19.7	19.8	20.4
PBF (%) <small>Percent Body Fat</small>	41.3	40.7	39.2	39.0	39.4	38.6	37.8	35.2
<input checked="" type="checkbox"/> Recent <input type="checkbox"/> Total	01.01.24 09:15	01.30.24 09:40	02.20.24 09:35	03.15.24 11:10	04.12.24 08:33	04.28.24 15:10	05.15.24 08:40	05.30.24 11:16

8 InBody Score

69 / 100 Points

* Total score that reflects the evaluation of body composition. A muscular person may score over 100 points.

9 Whole Body Phase Angle

$\phi(^{\circ})$ 50 kHz	4.3	4.4	4.2	4.1	4.3
	03.15.24 11:10	04.12.24 08:33	04.28.24 15:10	05.15.24 08:40	05.30.24 11:16

10 Weight Control

Target Weight	53.0 kg
Weight Control	- 6.1 kg
Fat Control	- 8.6 kg
Muscle Control	+ 2.5 kg

11 Waist-Hip Ratio

0.96	0.75 0.85
------	-----------

12 Visceral Fat Level

11	Low 10 High
----	-------------

13 Research Parameters

Fat Free Mass	38.3 kg (36.7 ~ 44.8)
Basal Metabolic Rate	1197 kcal (1255 ~ 1451)
Obesity Degree	112 % (90 ~ 110)
SMI	5.8 kg/m ²
Recommended calorie intake per day	1397 kcal


Calorie Expenditure of Exercise

Golf	104	Gateball	112
Walking	118	Yoga	118
Badminton	134	Table Tennis	134
Tennis	177	Bicycling	177
Boxing	177	Racketball	177
Mountain Climbing	193	Jumping Rope	207
Aerobics	207	Jogging	207
Soccer	207	Swimming	207
Japanese Fencing	295	Racketball	295
Squash	295	Taekwondo	295

* Based on your current weight

* Based on 30 minute duration

14 Impedance



Result Sheet Interpretation

1 Body Composition Analysis

Body weight is the sum of Total Body Water, Protein, Minerals, and Body Fat Mass. Maintain a balanced body composition to stay healthy.

2 Muscle-Fat Analysis

Compare the bar lengths of Skeletal Muscle Mass and Body Fat Mass. The longer the Skeletal Muscle Mass bar is compared to the Body Fat Mass bar, the better health status you are in.

3 Obesity Analysis

BMI is an index used to determine obesity by using height and weight. PBF is the percentage of body fat compared to body weight.

4 Segmental Lean Analysis

Analyze to evaluate if the muscles are adequately developed in the body segments. The top bar shows your muscle mass status compared to the ideal weight and the bottom bar shows your muscle mass status compared to your current weight.

5 Segmental Fat Analysis

Evaluate whether the amount of fat is adequately distributed in segments of the body. Each bar shows fat mass in comparison to the ideal amount.

6 Body Composition History

Measure your body composition periodically to monitor your progress.

7 Logo Customization

The Customized Logo can be applied on the Result Sheet. URL can also be applied at the bottom of the Result Sheet as well.

8 InBody Score

InBody Score is your score after evaluating your body composition.

9 Whole Body Phase Angle

Phase Angle is related to the health status of the cell membrane. Strengthening of the cellular membrane and structural function will increase the Phase Angle. In contrast, impairments to the cellular membrane can result in decreased Phase Angle.

10 Weight Control

Use the Target Weight, Weight Control, Fat Control, Muscle Control to set your own goal.

11 Waist-Hip Ratio (WHR)

Waist-Hip Ratio is the ratio of waist circumference to hip circumference.

12 Visceral Fat Level

Visceral Fat Level is an indicator based on the estimated amount of fat surrounding internal organs in the abdomen.

13 Research Parameters

Various research parameters such as Basal Metabolic Rate, Waist-Hip Ratio, Obesity Degree, Skeletal Muscle Mass Index (SMI), Body Cell Mass, and more are provided.

14 Impedance

Impedance is the resistance value measured when electrical currents are applied to the body. Based on the measured data, key body composition outputs can be analyzed. Impedance is also used for many research purposes.

*Additional InBody Results Sheet

InBody Result Sheet for Children, Thermal Result Sheet



1

InBody

[InBody270S]

InBody

inbody.com

ID	Height	Age	Gender	Test Date / Time
Jane Doe C	139.4cm	10	Male	06.21.2024 16:40

Body Composition Analysis

Total amount of water in my body	Total Body Water	(L)	19.2 (18.0 ~ 22.0)
What I need to build muscles	Protein	(kg)	5.1 (4.9 ~ 5.9)
What I need for strong bones	Minerals	(kg)	1.91 (1.66 ~ 2.04)
Where my excess energy is stored	Body Fat Mass	(kg)	8.8 (3.8 ~ 7.7)
Sum of the above	Weight	(kg)	35.0 (27.3 ~ 36.9)

Muscle-Fat Analysis

	Under	Normal	Over
Weight (kg)	55 70 85 100 115 130 145 160 175 190 205 %	35.0	
SMM (kg)	70 80 90 100 110 120 130 140 150 160 170 %	13.4	
Body Fat Mass (kg)	40 60 80 100 160 220 280 340 400 460 520 %	8.8	

Obesity Analysis

	Under	Normal	Over
BMI (kg/m ²)	7.9 10.9 13.9 16.4 18.6 20.2 22.2 24.2 26.2 28.2 30.2	18.0	
PBF (%)	0.0 5.0 10.0 15.0 20.0 25.0 30.0 35.0 40.0 45.0 50.0	25.2	

Growth Graph

Height : 50 ~ 75% Weight : 50 ~ 75%

Body Composition History

Height (cm)	136.5	137.2	138.6	139.4
Weight (kg)	35.1	35.6	37.3	35.0
SMM (kg)	13.3	13.0	12.9	13.4
PBF (%)	26.2	26.5	26.0	25.2
Recent <input checked="" type="checkbox"/> Total	09.10.23 09:15	11.30.23 09:40	01.02.24 09:35	06.21.24 16:40

Growth Score

86 / 100 Points

* If tall and within great body comparison standards, the growth score may surpass 100 points.

Obesity Evaluation

BMI	<input checked="" type="checkbox"/> Normal	<input type="checkbox"/> Under	<input type="checkbox"/> Slightly Over
PBF	<input type="checkbox"/> Normal	<input type="checkbox"/> Slightly Over	<input checked="" type="checkbox"/> Over

Research Parameters

Basal Metabolic Rate	936 kcal (948 ~ 1077)
Child Obesity Degree	109 % (90 ~ 110)

Results Interpretation

Growth Graph

Compares the height and weight among peers of the same age group.

Results Interpretation QR Code

Scan the QR Code to see results interpretation in more detail.

Impedance

Copyright © 1996~by InBody Co.,Ltd. All rights reserved. IR-English-UB-Child-A-240701

2

InBody

05/30/2024 11:16

ID : Jane Doe

Height : 156.9cm

Gender : Female

Age : 51

Weight : 59.1kg

Weight

59.1 kg

Muscle Mass

20.4 kg

It is the Mass of muscle attached to your bones. (Skeletal Muscle Mass)

Body Fat Mass

20.8 kg

Percent Body Fat

35.2 %

Normal Range (18.0~28.0)

Body Mass Index

24.0 kg/m²

Normal Range (18.5~25.0)

Basal Metabolic Rate

1197 kcal

it is the Minimum number of calories needed to sustain life at a resting state.

Waist Hip Ratio

0.96

Normal Range (0.75~0.85)

Visceral Fat Level

11

Normal Range (1~9)

Segmental Lean Analysis

1.92 kg

95.9 %

Normal

17.7 kg

97.4 %

Normal

2.00 kg

100.0 %

Normal

Left ————— Right

5.13 kg

81.0 %

Under

5.20 kg

82.2 %

Under

Segmental Fat Analysis

1.5 kg

168.6 %

Over

11.1 kg

223.3 %

Over

1.4 kg

163.1 %

Over

Left ————— Right

2.8 kg

124.1 %

Normal

2.8 kg

125.5 %

Normal

* Segmental fat is estimated.

InBody Score

69

Fat Control

- 8.6 kg

Muscle Control

+ 2.5 kg

Whole Body Phase Angle


4.3 °

Impedance

[000/000/000]

InBody Data Integration Solution

Manage and utilize your InBody data in various settings.



InBody Data Comprehension

Provide a health report to monitor your customers body composition goal.

Analytical Dashboard and Report

Get an intuitive analysis of your InBody data on the dashboard and see how your facility is operating with InBody.

Monitor Lifestyle Habits

Integrate InBody devices to monitor lifestyle habits and provide remote health management.

Access InBody Results Anywhere, Anytime

Through PC, tablet and smartphones, access your customer's InBody results anywhere, anytime.

API Integration


Upon customer consent, utilize InBody data through API and SDK.

Various File Formats

Print InBody data as an image, excel file etc.

Specifications

InBody270S Body Composition Analyzer



Bioelectrical Impedance Analysis (BIA) Measurement Items	Impedance (Z)	15 Impedance measurements by using 3 different frequencies (20 kHz, 50 kHz, 100 kHz) at each 5 segments of the body (Right Arm, Left Arm, Trunk, Right Leg, Left Leg)	Age Range	3+ years
	Phase Angle (Ø)	1 Phase Angle Measurements by Using 1 Frequency (50 kHz) at Whole Body	Height Range	95 - 220 cm (3 ft 1.4 in - 7 ft 2.6 in)
Electrode Method	4 electric poles 8 Points Touch type electrode measurement		Language Support	InBody supports over 30 languages.
Measurement Method	Direct Segmental Multi-frequency Bioelectrical Impedance Analysis Method, DSM-BIA Method		Outputs (InBody Result Sheet)	Results and Interpretations <ul style="list-style-type: none">· Body Composition Analysis (Total Body Water, Protein, Minerals, Body Fat Mass, Weight)· Muscle-Fat Analysis (Weight, Skeletal Muscle Mass, Body Fat Mass)· Obesity Analysis (Body Mass Index, Percent Body Fat)· Segmental Lean Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg)· Segmental Fat Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg)· Body Composition History (Weight, Skeletal Muscle Mass, Percent Body Fat)· InBody Score· Whole Body Phase Angle (History)· SMI (History)· Weight Control (Target Weight, Weight Control, Fat Control, Muscle Control)· Nutrition Evaluation (Protein, Minerals, Fat Mass)· Obesity Evaluation (BMI, Percent Body Fat)· Body Balance Evaluation (Upper, Lower, Upper-Lower)· Waist-Hip Ratio (Graph)· Visceral Fat Level (Graph)· Research Parameters (Skeletal Muscle Mass, Fat Free Mass, Basal Metabolic Rate, Waist-Hip Ratio, Waist Circumference, Visceral Fat Level, Obesity Degree, FFMI, FMI, SMI, SMM/WT)· Recommended calorie intake per day· Calorie Expenditure of Exercise· Sarcopenia Parameter (SMI, HGS)· Blood Pressure (Systolic, Diastolic, Pulse, Mean Artery Pressure, Pulse Pressure, Rate Pressure Product)· QR Code· Results Interpretation QR Code· Whole Body Phase Angle (50 kHz)· Impedance (Each segment and each frequency)
Body Composition Calculation	No use of Empirical Estimation		Outputs (InBody Result Sheet for Children)	Results and Interpretations <ul style="list-style-type: none">· Body Composition Analysis (Total Body Water, Protein, Minerals, Body Fat Mass, Weight)· Muscle-Fat Analysis (Weight, Skeletal Muscle Mass, Body Fat Mass)· Obesity Analysis (Body Mass Index, Percent Body Fat)· Growth Graph (Height, Weight, BMI)· Body Composition History (Height, Weight, Skeletal Muscle Mass, Percent Body Fat)· Whole Body Phase Angle (History)· SMI (History)· Growth Score· Weight Control (Target Weight, Weight Control, Fat Control, Muscle Control)· Obesity Evaluation (BMI, Percent Body Fat)· Nutrition Evaluation (Protein, Minerals, Fat Mass)· Body Balance (Upper, Lower, Upper-Lower)· Research Parameters (Skeletal Muscle Mass, Fat Free Mass, Basal Metabolic Rate, Child Obesity Degree, FFMI, FMI, SMI, SMM/WT)· Sarcopenia Parameter· Blood Pressure (Systolic, Diastolic, Pulse, Mean Artery Pressure, Pulse Pressure, Rate Pressure Product)· QR Code· Results Interpretation QR Code· Whole Body Phase Angle (50 kHz)· Impedance (Each segment and each frequency)
Logo Display	Name, Address, and Contact Information can be shown on the InBody Results Sheet.		Outputs (InBody Thermal Result Sheet)	Total Body Water, Protein, Minerals, Weight, Muscle Mass, Body Fat Mass, Percent Body Fat, BMI, Basal Metabolic Rate, Waist-Hip Ratio, Waist Circumference, Visceral Fat Level, FFMI, FMI, SMI, SMM/WT, Segmental Lean Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg), Segmental Fat Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg), InBody Score, Fat Control, Muscle Control, Whole Body Phase Angle (50 kHz), Impedance (Each segment and each frequency)
Digital Results	LCD Monitor, Data management software LookinBody120			
Type of Results Sheets	InBody Results Sheet, InBody Results Sheet for Children, InBody Thermal Results Sheet			
Data Storage	Test results can be saved if the member ID is utilized. The InBody can save up to 100,000 results.			
Test Mode	Self Mode, Professional Mode			
Administrator Menu	<ul style="list-style-type: none">· Setup: Configure settings and manage data· Troubleshooting: Additional information to help use the InBody			
USB Thumb Drive	Copy, Backup or restore the InBody test data (data can be viewed on Excel or LookinBody data management software)			
Barcode Reader	The member ID will be automatically inputted when the barcode ID is scanned.			
Backup data	Backup data saved in the InBody by using a USB Thumb Drive, Restore results on the InBody from a backup file.			
Applied Rating Current	300 µA (± 30 µA)			
Adapter ② (DELTA)	Power Input	AC 100 - 240 V, 50 - 60 Hz, 1.5 A - 0.75 A		
	Power Output	DC 12 V, 5.0 A		
Adapter ① (MEAN WELL)	Power Input	AC 100 - 240 V, 50 / 60 Hz, 1.0 A - 0.5 A		
	Power Output	DC 12 V, 3.34 A		
Display Type	480 × 800 7 inch Color TFT LCD			
Internal Interface	Touchscreen, Keypad			
External Interface	RS-232C 1EA, USB HOST 2EA, USB SLAVE 1EA, LAN (10/100T) 1EA, Bluetooth 1EA, Wi-Fi(2.4G/5G) 1EA			
Compatible Printer	Laser/Inkjet PCL 3 or above and SPL			
Dimension	356 (W) X 796 (L) X 995 (H) : mm 14.0 (W) X 31.3 (L) X 39.2 (H) : inch			
Device Weight	13.4 kg (29.5 lb)			
Testing Duration	About 30 sec.			
Operation Environment	10 - 40 °C (50 - 104 °F), 30 - 75 % RH, 70 - 106 kPa			
Storage Environment	-10 - 70 °C (14 - 158 °F), 10 - 80 % RH, 50 - 106 kPa (No Condensation)			
Testing Weight Range	2 - 250 kg (4.4 - 551.2 lbs)			

* Specifications may change without prior notice.
* "QR Code" is registered trademark of DENSO WAVE INCORPORATED.

InBody

InBody HQ [KOREA]
InBody Co., Ltd.
InBody Bldg., 625, Eonju-ro, Gangnam-gu, Seoul
06106 Republic of Korea
TEL: +82-2-501-3939 FAX: +82-2-6919-2417
Website: inbody.com
E-mail: info@inbody.com

InBody China [CHINA]
Biospace China Co., Ltd.
903/904, XingDi Plaza, No.1698 YiShanRoad,
Shanghai 201103 China
TEL: +86-21-6443-9705 FAX: +86-21-6443-9706
Website: inbodychina.com
E-mail: info@inbodychina.com

InBody Oceania [AUSTRALIA]
InBody Oceania Pty Ltd.
U2/82-86 Minnie Street, Southport, Queensland
TEL: +61-7-5681-1900
Website: au.inbody.com
Email: oceania@inbody.com

InBody USA [USA]
Biospace Inc. dba InBody
13850 Cerritos Corporate Dr. Unit C Cerritos,
CA 90703 USA
TEL: +1-323-932-6503 FAX: +1-323-952-5009
Website: inbodyusa.com
E-mail: info.us@inbody.com

InBody Europe B.V. [NETHERLANDS]
InBody Europe B.V.
Gyrocoopweg 122, 1042 AZ, Amsterdam,
The Netherlands
TEL: +31-20-238-6080 FAX: +31-6-5734-1858
Website: nl.inbody.com
E-mail: info.eu@inbody.com

Biospace Latin America [MEXICO]
Biospace Latin America S. de R.L. de C.V.
Insurgentes Sur 1457, Piso 15 Int.2. Col. Insurgentes
Mixcoac, Alcaldia Benito Juarez, C.P. 03920, Ciudad
de México, México
TEL: +52-55-5025-0147
Website: inbodymexico.com
E-mail: info.mx@inbody.com

InBody BWA Inc. [USA]
InBody BWA Inc.
2550 Eisenhower Avenue, Suite C 209, Audubon, PA 19403
TEL: +1-610-348-7745
Website: inbodybwa.com
E-mail: bwainquiries@inbody.com

InBody Europe B.V. Niederlassung Deutschland [GERMANY]
InBody Europe B.V.
Mergenthalerallee 15-21, 65760 Eschborn, Germany
TEL: +49-6196-76-916-62 FAX: +49-6196-76-916-11
Website: de.inbody.com
E-mail: erfolg@inbody.com

InBody Asia [MALAYSIA & SINGAPORE]
InBody Asia Sdn. Bhd.
Unit 3A-11, Oval Damansara, 685 Jalan Damansara
Kuala Lumpur, WP KL 60000 Malaysia
TEL: +60-3-7732-0790 FAX: +60-3-7733-0790
Website: inbodyasia.com
E-mail: info@inbodyasia.com

InBody Japan [JAPAN]
InBody Japan Inc.
Tani Bldg., 1-28-6, Kameido, Koto-ku, Tokyo
136-0071 Japan
TEL: +81-3-5875-5780 FAX: +81-3-5875-5781
Website: inbody.co.jp
E-mail: inbody@inbody.co.jp

InBody UK [UNITED KINGDOM]
11 Phoenix Park, Telford Way, Stephenson Industrial
Estate, Coalville LE67 3HB, United Kingdom
TEL: +44-1530-569620
Website: uk.inbody.com
E-mail: uk@inbody.com

InBody India [INDIA]
InBody India Pvt.Ltd.
57/57 A, 1st floor, Raj Industrial Complex, Military Road,
Marol, Andheri (East), Mumbai- 400059, Maharashtra,
India
TEL: +91-22-6223-1911
Website: inbody.in
E-mail: india@inbody.com



For more details about the patents that we acquired, please visit our website or refer to the patent gazette of intellectual property office of each country. (Korea, U.S, China, Japan)