


InBody770s

INBODY HAS THE ABILITY TO TELL YOU WHAT YOUR BODY CANNOT TELL, IN 30 SECONDS!

ESSENTIAL PARAMETERS MEASURED

1. VISCERAL FAT AREA/LEVEL

- Body fat within the abdominal cavity
- Main cause of diabetes & heart attack
- Normal: $< 100 \text{ cm}^2$ – Healthy
- Over: $100 \sim 150 \text{ cm}^2$ – Moderate risk for Type 2 Diabetes & Heart Disease
- Extremely over: $> 150 \text{ cm}^2$ – High risk for Type 2 Diabetes & Heart Disease




2. PERCENTAGE BODY FAT

- Body fat mass relative to total body weight
- Unhealthy zone:
Females: **23%** and above
Males: **18%** and above



3. WAIST HIP RATIO (WHR)

- Waist circumference relative to hip circumference
- Ideal WHR:
Females: **0.8**
Males: **0.9**



5. SKELETAL MUSCLE MASS

- Weight of muscles attached to your skeleton
- Skeletal muscle health and mass varies in individuals
- Should be in the normal or over range



7. ECW RATIO ANALYSIS

- Indicates whether body water is balanced
- **0.36 ~ 0.39** – healthy
- **< 0.36** – dehydrated
- **> 0.39** – edema (water retention)



9. BASAL METABOLIC RATE (BMR)


- Amount of calories burnt by the body at rest
- Females: **1200kcal**
- Males: **1500kcal**



4. WHOLE BODY CELL PHASE

ANGLE @ 50 kHz

- How well protected are your cells
- $> 4^\circ$: normal
- $< 4^\circ$: unhealthy individual body immunity level
- $< 3^\circ$: critically unhealthy




6. WEIGHT CONTROL

- How much muscles and fats to be gained or lost
- Weight control varies in individuals



8. BONE MINERAL CONTENT

- How rich your bones are in minerals
- Bone mineral content can indicate **osteoporosis**



10. INBODY SCORE

- Total score that reflects your body's physical health
- The **more muscles** you have, the **higher** your InBody score!



Segmental Lean Analysis

	Under	Normal	Over
Right Arm (kg)	2.55	2.96	3.37
Left Arm (kg)	2.55	2.96	3.37
Trunk (kg)	21.85	23.9	25.95
Right Leg (kg)	6.79	7.79	8.79
Left Leg (kg)	6.79	7.79	8.79

Ideal lean value: 100% or more

11. SEGMENTAL LEAN ANALYSIS & SEGMENTAL FAT ANALYSIS

- Show distribution of **muscle mass** & **fat** between arms, legs and trunk
- Imbalance of arms and legs lean and fat mass is not ideal
- The result varies among individuals
- Refer to the scale given under SMM(kg) & Body Fat Mass (kg) in % for the normal range

Segmental Fat Analysis

Right Arm	(0.5 kg)	77.1%
Left Arm	(0.4 kg)	70.9%
Trunk	(6.2 kg)	149.8%
Right Leg	(1.8 kg)	107.0%
Left Leg	(1.8 kg)	106.9%

Ideal fat value: 100% or less

TO KNOW MORE ABOUT YOUR BODY, TAKE THE INBODY TEST NOW!



Health is Wealth

DISTRIBUTOR:
SMITECH (ASIA) PTE LTD

75 BUKIT TIMAH ROAD, #06-13 BOON SIEW BUILDING
SINGAPORE 229833 TEL: +65 6339 5277 FAX: +65 6339 0253

E: sales@smitechasia.com.sg or smitech@pacific.net.sg

W: www.smitechasia.com.sg



Your Partner
For Healthy Lifestyle