

Oxygen Concentrator OLV-10



Two colors available



Standard Features:

- France imported CECA molecular sieve;
- Long lifespan more than 20000 hours;
- 2 years warranty;
- CE ISO approved;



Standard Accessories

NO.	Name	Q'ty(pcs)	NO.	Name	Q'ty(pcs)
2-1	Nasal Pipe	2	2-4	Humidifier Bottle	1
2-2	Dust Filter	1	2-5	Power Cord (EU/US/CN/AU/UK)	1
2-3	HEPA Filter	1	2-6	User Manual	1

Standard Parameters

Model NO.	OLV-10	OLV--5
Flow Rate	1-10L/MIN	1-5L/MIN
Concentration	93% ± 3 (V/V)	93% ± 3(V/V)
Power Consumption	610W	<u>330W</u>
Sound Level	≤50dB	<u><45dB</u>
Outlet Pressure	40-80 kPa	<u>40-80 kPa</u>
Net Weight	24KG	<u>15KG</u>
Dimensions	372×340×612 mm	<u>372x340x612 mm</u>
Carton Size	445×400×660 mm	<u>445x400x660 mm</u>

Option: SMITECH Oxygen Concentrator OLV-5 5 Litre/Min Capacity of Pure Oxygen for 10,000 hours
 SMITECH Oxygen Warmer with temperature control, thermometer and 2 oxygen masks
 SMITECH Handheld Pulse Oximeter with SPO2 and Pulse Rate LCD Display

Exclusive Distributor: **SMITECH (ASIA) PTE LTD**

75 BUKIT TIMAH ROAD #06-13 BOON SIEW BUILDING
 SINGAPORE 229833

Tel: +65 6339 5277 Fax: +65 6339 0253 Web: www.smitechasia.com.sg

Email: sales@smitechasia.com.sg; smitech@pacific.net.sg;